## ENGLISH

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

**IF YOU ARE USING THE CRYO/CUFF WITH THE AIRCAST CRYO/CUFF IC UNIT BE SURE YOU READ AND UNDERSTAND ALL INSTRUCTIONS LOCATED IN THE AIRCAST CRYO/CUFF IC INSTRUCTIONS FOR USE.**

**INTENDED USER PROFILE:** The intended user should be a licensed medical professional, the patient, the patient's caregiver or a family member providing assistance. Users should be able to:  
• Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

**OPERATING PRINCIPLE:** The AirCast Cryo/Cuff provides cold therapy and compression.

**INTENDED USE/INDICATIONS:** The Aircast Cryo/Cuff combines focal compression with cold to minimize swelling and pain.

**FREQUENCY OF USE:** The length of use and frequency of use of the Cryo/Cuff are determined by the treating physician and/or individual's behavior and the patient. The user should be advised to use serious injury be sure to check skin frequently.

**CONTRAINDICATIONS:** Cryotherapy should not be used on persons with Raynaud's or other vasospastic disease, cold hypersensitivity, decreased skin sensitivity, compromised local circulation, diabetes, stroke, cerebral aneurysm, arterial hypertension, hyperostosis osteoarthritis, osteoarthritis, distal peripheral vascular disease causing ischemia or poor local circulation, local tissue infection, cold allergy, cold hypersensitivity, previous cold injury and paronychia cold hemoglobinuria.

Do not use if you are allergic to any of the materials contained within this product.

### WARNINGS AND PRECAUTIONS:

This device can be cold enough to cause serious injury, including tissue damage. You must be able to check your skin condition under the device and if you cannot check your skin condition frequently. People are sensitive to cold in diverse ways and may react differently to cold treatment. Check for increased pain, burning, welts, numbness, tingling, increased redness, discoloration, itching, increased swelling, blisters, irritation or other changes in skin condition under the Cuff or around the treatment area. If you experience any of these conditions, immediately discontinue use of the device and contact your physician.

When using the device, you must occur every 1-2 hours on a routine basis. Caution must be taken by routinely monitoring the patient during use on children, elderly, incapacitated patients, severe cardiovascular diseases, those with decreased skin sensitivity or circulation due to poor nutrition, malnutrition, dehydration, hypotension, hypoxemia, narcotic, tobacco and alcohol use. If such monitoring is to occur outside the presence of the licensed healthcare professional, he or she must ensure that the patient or patient caretaker understands how such monitoring is to be conducted.

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When applying the device wrap to the skin, apply a barrier (ie tape or cloth) between the wrap and your skin, and do not let the wrap sit around the armpits. No part of the wrap should touch the skin for any period of time.

Check for moisture. If moisture is present on the Cuff, immediately discontinue use of the device. Do not use an elastic wrap in conjunction with the Cryo/Cuff.

Dressings used under the Cryo/Cuff should be applied lightly.

When filling the Cryo/Cuff, do not raise the cooler higher than 15 inches (38 cm) above the cuff to avoid excessive pressure.

Reduce pressure with any sense of discomfort, numbness or tingling of the limb.

The Cryo/Cuff is designed for single patient use only and may be used on the same patient for the length of treatment as permitted.

The device should be cleaned if it is used on the same patient for an extended period of time (see "Cleaning Instructions" for cleaning information).

Keep the device and packaging materials away from small children and animals. These items pose a risk for suffocation or strangulation.

It could be unsafe to use accessories, detachable parts and materials, or interconnect to other accessories unless described in these instructions, or otherwise modified equipment.

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**PREPARE COOLER** - Connect a built tube to cooler. Add water to the inside cooler, then add ice to top of cooler. Lay the insulation cloth on top of the ice. Attach cooler lid snugly. Allow 5 minutes with occasional shaking to chill water.

**APPLY THE EMPTY CUFF** - Center cuff over injured area. Aircast logo must be upright facing out. Secure strap (Fig. 1).

**FILL CUFF** - Connect blue tube to the cuff. Insert for a click/snag. When using the Gravity Cryo/Cuff (non-motorized), open air vent on cooler. Insert a minor or two for the water to mix with the ice and reach the cuff. Then close the cooler and repeat the filling process (see step 3).

**RECHILL WATER** - Reconnect the blue tube to cuff, lower the cooler below cuff and the warmed water will drain from the cuff into the cooler. Allow a minute or two for the water to mix with the ice and reach the cuff, then close the cooler and repeat the filling process (see step 3).

**CLEANING INSTRUCTIONS:** After use, completely drain water from cuff, tube, and cooler. (Do drain tube, elevate tube while pressing tip of quick-disconnect.) Store cooler with top off to allow drying. Periodically clean cuff, tube and cooler using a few ounces of liquid soap added to hot water cooler. Recycle soap/water into a drain between cooler and cuff a few times, then repeat with clean water only. Rinse, clean, dry, and store.

**MATERIAL CONTACTS:** Polyurethane/Nylon, Nylon, Polyether foam.

**WARRANTY:** DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

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**NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.**

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