DONJOY®

Support Everest™ Ice II

Support Everest™ Ice II

Hinged knee brace with Cold Compression

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

SUPPORT EVEREST™ II

SUPPORT EVEREST™ ICE II

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

SUPPORT EVEREST™ II HINGED KNEE BRACE

INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to:

- Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS

- Post-traumatic, pre-op or post-operative use, functional rehabilitation use, return to activities.
- Light to moderate ACL/PCL injuries, moderate to severe lateral ligament injuries (MCL, LCL), meniscal injuries, knee surgery.

CONTRAINDICATIONS

None

APPLICATION INSTRUCTIONS

The knee brace has 2 lateral hinges blocked in extension. It also has a removable posterior stay along the popliteal area. This stay is pre-shaped, but can be moulded according to patient's needs or it can be removed.

With the patient lying on their back, open the brace.

- 1 - Position the patient's leg on the brace, with the knee centred between the 2 bilateral hinges. Close the anti-slip strap on the calf using the hook and loop closure.
- 2 - Secure the brace by closing the thigh cuff and calf cuff using the hook and loop closures.
- 3 - Adjust the 4 straps to the required length and close using the buckle/clip closure. Tighten the brace according to level of compression required.

Hinges

- To adjust flexion/extension, open the cap of each hinge, position the flexion/extension stops at the desired angle, and close the hinge cap.
- To avoid damaging the hinges, the flexion/extension stops must be positioned at the same angle on both hinges.

WARNING & PRECAUTIONS

- We recommend that the first application be made under supervision of a health care professional.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.

CARE

- Hand wash in lukewarm water (30 ° C) with a mild detergent
- Air dry (do not place on radiators/heaters or other heat sources)
- Do not use bleach or chlorine, do not iron.

KNEE CUFF

- Wrap may be hand washed and air dried. Gel inserts may be cleaned with water, moist towel or sponge.

COMPOSITION

Knee brace: Aluminum, polyester, polyamide, monofilament, polyurethane, cotton, elastomer, resin, acetal (POM), polypropylene, polyethylene

Knee Cuff –

Inflatable Pad: Polyester, nylon

Pump: Silicone, Polyethylene, PVC, ABS Resin (thermoplastic polymer), synthetic rubber and steel wire

Cold Pack - Gel: propylene glycol, sodium carboxymethyl cellulose, water.

Fabric: cotton, nylon

CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS.

INTENDED FOR USE ON A SINGLE PATIENT

WARRANTY DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

SUPPORT EVEREST™ ICE II HINGED KNEE BRACE WITH COLD COMPRESSION

For treatment combining Cold therapy and Compression, the Knee Brace is delivered with a KNEE CUFF containing 1 Inflatable Pad, 1 Hand Pump and 2 Cold Packs.

INTENDED USER PROFILE

The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to:

- Read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE/INDICATIONS

- To be used ONLY in conjunction with SUPPORT EVEREST Knee Brace.
- Localized external application of cold therapy and compression after surgery or injury, to help minimize inflammation, swelling, and pain.
- Providing immobilization or controlled movement of the limb or body segment.

CONTRAINDICATIONS

DO NOT use this device on patients with:

Raynaud's or other vasospastic diseases; diabetics; sickle-cell anemia, arteriosclerosis, cryoglobulinemia, hypercoagulable clotting disorders or other peripheral vascular diseases causing ischemia or poor local circulation; compromised local circulation; local tissue infection; cold allergy, cold hypersensitivity, previous cold injury.

APPLICATION INSTRUCTIONS

Prior to use, detach Cold Pack from Knee Cuff and place in freezer for 4 hours minimum. (Domestic deep freeze temperature between -12°C & -20°C).

A - Once Cold Pack is cold / frozen, attach the Cold Pack onto the Knee Cuff, using the hook and loop.

Feel Cold Pack and ensure temperature is adequate before use.

• **WARNING!** Be sure to place a barrier/cloth between your skin and the Cold Pack. (See "WARNING & PRECAUTIONS"). Cold Pack temperature should not be lower than 5C.

B - Place the Knee Cuff around the knee, with the patella opening centred on the knee cap. Secure the Knee Cuff by attaching the Thigh Strap and the Calf Strap using the hook and loop closures.

C - Position the SUPPORT EVEREST KNEE BRACE above the Knee Cuff. (See "APPLICATION INSTRUCTIONS" 1 to 3).

D - Screw the pump to the Knee Cuff, turn the valve to 'ON' position. Inflate the Air Cell to the desired level of compression.

Once the desired inflation level has been reached turn the valve to 'OFF' and remove the pump.

• **WARNING!** Do not inflate and tighten the Knee Cuff excessively. Ensure level of compression is not too tight so as not to restrict blood or fluid flow. (See "WARNING & PRECAUTIONS").

• **WARNING!** Upon instruction of a Health Care Professional: Apply the Knee Cuff to injured area for time recommended by a Health Care Professional. Do not exceed 20-30 minutes period of application. Alternate with adequate rest periods.

If needed, renew application using second Cold Pack. Replace first Cold Pack in freezer until needed. Do not sleep with the Knee Cuff on.

E - To deflate, screw the pump to the Knee Cuff, position the valve on ON, and press the button at the base of the pump to release the air.

WARNING & PRECAUTIONS

WARNING! THIS DEVICE CAN BE COLD ENOUGH TO CAUSE SERIOUS INJURY. SERIOUS ADVERSE REACTIONS AND SAFETY HAZARDS MAY OCCUR WHEN USING THIS DEVICE.

USE OF KNEE CUFF WITH INFLATABLE PAD & COLD PACK :

- A health care practitioner must state how long and how often the device should be used and the length of breaks between uses.
- Extreme care must be taken when using any cold therapy system. Any cold products may cause injuries, e.g. burns, blisters, frostbite, swelling, discoloration and skin or tissue necrosis if improperly used or monitored.
- **WARNING!** DO NOT inflate and tighten the Knee Cuff excessively. Ensure level of compression is not too tight so as not to restrict blood or fluid flow.
- **WARNING!** This device can be cold enough to cause serious injury, including tissue necrosis. You must be able to check your skin condition under the Cold Pack frequently (at least every hour). DO NOT use this device if you cannot check your skin condition frequently (at least every hour). If you develop an allergic reaction and/or experience itchy, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately. Check for increased pain, burning, numbness, tingling, increased redness, discoloration, itching, increased swelling, blisters, irritation or other changes in skin condition under the Cold Pack or around the treatment area. If you experience any of these conditions, immediately discontinue use of this device and contact your physician.
- **WARNING!** Application of the Cold Packs directly on the skin may result in serious injury, including tissue necrosis. DO NOT let any part of the Cold Pack touch your skin. ALWAYS use with a barrier/cloth between your skin and the Cold Pack.
- Keep the product away from children.
- All Cold Pack ingredients are non-toxic. They are, however, not intended for ingestion.
- Do not heat the Cold Packs (microwave, oven, hot water, etc...)
- Do not use on open wound(s) and/or sensitive skin.
- Do not use this product if it is damaged or defective.
- Do not use this product in any way that may damage its integrity.
- Do not sit or lay over the Cold Pack or fall asleep with the Cold Pack on the skin.
- If contents get in contact with eye or skin, rinse with water immediately.

LIMITATIONS ON SPECIAL PATIENT POPULATIONS

Limit the use of this device with patients who are unresponsive, incapacitated, have altered mental status or altered pain perception. Post-surgical patients under sedation or on analgesics or anesthetics, as well as patients taking hypnotics, anxiolytics, or antidepressants, must be monitored frequently during use of this device. These patients may not be able to perceive pain, burning, numbness, tingling or decreased sensation and may be susceptible to injury. Discontinue cold therapy immediately at the first sign of cold injury.

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MDSS GmbH
DJO FRANCE S.A.S.