ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a licensed medical professional, the patient, the patient's caretaker, or a family member providing assistance. The user should be able to read, understand and be physically capable to perform all the directions, warnings and cautions provided in the information for use.

INTENDED USE / INDICATIONS: The Air-Stirrup provides support to prevent further injury by stabilizing the ankle.

PERFORMANCE CHARACTERISTICS: Soft-good/semi-rigid designed to restrict motion through elastic or semi-rigid elastication.

CONTRAINDICATIONS: N/A

WARNINGS AND PRECAUTIONS:

Like all lower extremity devices, such as casts or braces, patients without sensation (i.e. post-op anesthesia, neuropathies, etc.) should be monitored frequently for "hot spots", skin irritation or wound management.

Always wear an absorbent sock or ankle wrap when using the Air-Stirrup Ankle Brace.

- Do not use over open wounds.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, you should contact your doctor immediately.
- Do not use this device if it was damaged and/or package has been opened.
- If you develop an allergic reaction and/or experience itching, red skin after coming into contact with any part of this device, please stop using it and contact your healthcare professional immediately. **NOTE:** Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INFORMATION:

1. ADJUST HEELPAD WIDTH

- Open brace and lift bottom of aircells until heel straps are fully exposed.
- Peel up heel straps and adjust heelpad width for a snug fit (Fig. 1).
- Firmly press heel straps and aircells back in place.

2. APPLY BRACE

- Put on absorbent sock or ankle wrap.
- Place round edge of heelpad under heel.
- Align brace sides with ankle (Fig. 2).

3. SECURE STRAPS

- Secure straps from bottom to top.
- Put on lace shoe or sneaker.

4. ADJUST BRACE

- Squeeze brace sides together with one hand.
- Tighten straps from bottom to top with other hand (Fig. 3).
- Tighten for comfortable support.

ANKLE BRACE FIT ADJUSTMENT

Ankle Brace aircells are preinflated and normally do not require adjustment.

If you experience pinching or uneven pressure:

- Reposition brace or,
- Readjust heelpad and straps. Use lines on heelpad as a reference.

If condition persists, adjust aircell volume.

(NOTE: Too much air will reduce support)

- Unfold valve on top of aircell.
- Insert inflation tube into valve 1/2 inches.
- Remove a small amount of air by gently squeezing aircell or.
- Add a small amount of air by blowing into tube.
- Squeeze valve flat just below tube, and slowly remove tube.
- Fold valve.

HIGH ALTITUDE:

At high altitudes the aircells will expand beyond their optimal level, and air will have to be removed from the aircells (see ankle brace fit adjustment). When flying, readjust straps to a comfortable pressure.

CARE:

The entire Ankle Brace can be hand washed in lukewarm water with mild soap and air dried.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

FOR SINGLE PATIENT USE ONLY.

~~NOT MADE WITH NATURAL RUBBER LATEX.~~

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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