ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE: The intended user should be a licensed medical professional, the patient or the patient's caregiver. The user should be able to read, understand and be physically capable of performing the directions, warnings and cautions in the information for use.

INDICATIONS: The Procure Nextep Contour W/Air walker is designed to provide support and immobilization to the lower leg, ankle and foot. It may be suitable for use following stress fractures of the lower limb, stable ankle fractures, Achilles tendon repair and acute lateral ankle sprains. Providing immobilization or controlled movement of a limb or body segment.

CONTRAINDICATIONS: This device is contraindicated for unstable fractures of the lower leg, ankle and foot. Do not use if you are allergic to any of the materials contained within this product.

WARNINGS AND PRECAUTIONS:

- This device is to be used under the supervision of a healthcare professional. The determination of when to apply the brace and the frequency and duration of use should be strictly at the discretion of the treating physician. Always consult with your physician before making changes to the brace.
- Application of this device is recommended only when the fracture is demonstrably stable and there are acceptable limits of angular and rotational deformity.
- Do not use this device on patient's incapable of communicating physical discomfort.
- Patients with diabetes, vascular deficiencies and neuropathies should be advised to perform circulation checks and consult your physician for further information.
- Accessory pads may be used where necessary to aid in patient comfort.
- Do not over-inflate aircells. Improper aircell inflation may cause significant skin irritation, especially in people with diabetic neuropathy. Reduce aircell pressure with any visual skin changes or reported discomfort.
- Do not over tighten straps. This may result in reduction of blood flow or sensation.
- Patient's without sensation (i.e. post-op anesthesia, neuropathies, etc.) should be monitored frequently for skin damage, skin irritation or wound management.
- Use caution when walking on slippery or wet surfaces to avoid injury.
- Do not use over open wounds.
- Do not modify the device or use the device other than intended.
- At high altitudes the aircells will expand beyond their optimal level. Adjust aircell support by tightening or loosening straps.
- Do not use this device if it is damaged and/or the packaging has been opened.
- If pain, swelling, changes in sensation or other unusual reactions occur while using this product, consult your medical professional immediately.

NOTE: Contact manufacturer and competent authority in case of a serious incident arising due to usage of this device.

APPLICATION INFORMATION:

Inflate and deflate bladder once before use to ensure integrity of air chambers.

1. Select correct size of walker for the patient. The end of the patient's toes should be within 1 in (2.5 cm) from, without extending past, the front end of the walker. Loosen straps and remove liner from boot. (Fig.1)
2. Place foot in liner and secure with contact closure. Make sure the heel fits snugly into the posterior portion of the liner. Secure the liner firmly around the foot, ankle, and lower leg. (Fig.2)
3. Lift the end of the boot inner sole and position one end of the toe wrap on the contact closure strip. (Fig.3)
4. Step into boot with heel seated firmly to the back of the shell. Fold the toe wrap over toes and secure on liner. (Fig.4)
5. Position the front shell so the distal edges (toe region) fit inside the rear shell and the proximal edges (shin region) fit outside the rear shell. Fasten the front shell to the back shell by securing the straps, starting with the instep strap first, followed by all other straps. Top strap must be secured outside of shell. (Fig.5)
6. To inflate walker liner press and release the air bulb several times until liner is providing comfortable cushion and support. DO NOT OVER INFLATE. (Fig.6a)
7. To release air from the walker liner, turn the inflation/deflation knob several times counter-clockwise. When re-inflating, turn knob clockwise until snug. Release air if liner feels too tight. (Fig.6b)

USE AND CARE: To clean liners, hand wash in cold water using mild detergent and air dry. Note: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

Always consult with your physician or therapist before making changes to the brace.

MATERIAL COMPONENTS:

- Shell: POLYPROPYLENE
- Outer Sole: NATURAL RUBBER
- Inner Sole: EVA
- Straps: NYLON
- Liner: NYLON AND FOAM

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale. To the extent the terms of this warranty are inconsistent with local regulations, the provisions of such local regulations will apply.

INTENDED FOR SINGLE PATIENT USE.

NOT MADE WITH NATURAL RUBBER LATEX.
RX only.

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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